


Hammarplast Medical

Disposables



HpmTM
For better care!

Hpm Products

We challenge ourselves to continuously improve. We are guided by our expertise, curiosity and creativity of finding new environmentally friendly medical solutions. In close co-operation with our customers, we want to contribute to the best possible medical care.

Class 1M medical measures - for your safety

We would like to inform you about some very important rules that exist for the medicine cup regarding how the product should be marked in order to fulfil the Medical Device regulations 93/42/EEC for dispensing cups used for both liquid medicine and tablets on the European market.

Our medicine cup 30 ml, used both for fluid medicine and tablets, is classified as a class 1M product.

Class 1M means that there is a measurement function with a scale on the product. The scale is designed to measure liquid medicine so that the correct volume is distributed to the patient.

All products intended to measure are classified as class 1M (The product has a measuring function). According to the European Medical Device regulation 93/42/EEC Annex VII, paragraph 5 of MDD, it is required that the manufacturer of class 1 devices with a measuring function must follow the regulations of the procedures referred to in annex IV, V and VI.



You can read in Annex VI:

“The manufacturer affixes the CE marking in accordance with article 17 and draws up a written declaration of conformity. The CE marking must be accompanied by the identification number of the Notified Body which performs the task referred to in this Annex.”

A product classified as class 1M product must have Notified Body mark and CE mark on the product.

CE
marked

Notified Body
marked

The Notified Body Number (N.B) signifies that the product is following the regulations of class 1M, where the Notified Body is inspecting the manufacturing company regarding the calibration of the scale.

Our Notified Body number 0413, implies that our Notified Body Intertek, inspects our calibration system of the medicine cups and controls that we follow the manufacturing rules of the product according to the Medical Device regulations.

An exact scale is very important when you distribute medicine to a patient. Note that the medicine description always prescribes an exact dose. If the product is not marked with Notified Body number the scale can not be trusted.

If the patient gets an overdose of medicine the risk of side effects increases, and if you give the patient too little medicine the desired effect reduces, which may pose discomfort to the patient.

This is especially important with small children where the given dose of medicine is based on the weight of the child to be exact.

Hammarplast Medical manufactures medical devices and is well educated in the Medical Device Regulations. We ensure our customers that the products provided by Hammarplast Medical are in order according to the MDD regulations.



Products



Medicine measures and lids



ECO+ a new range of products from Hpm™ made from renewable raw material.

The Hpm™ vision calls for us to support the transformation towards a more sustainable society. One of our main goals is to increase production and sales of products made of renewable materials in order to reduce our carbon footprint. Launching ECO+ medicine measure is one of our important contributions for a sustainable world.



ECO+ medicine measure

30 ml medicine measure graduated each ml.

The measures are made of polypropylene and the writable lids of white polyethylene.

Measure: staple/carton 90/9 000

Lid: bag/carton 500/9 000

Art no:

02350 nature

02351 yellow

02352 blue

02353 red

02354 green

00364 lid



ECO+ a new medicine measure from Hpm™ made of renewable raw material.

In our strategy for a more sustainable society, we have chosen to manufacture our products according to the mass balance system. Yearly we will be audited by the independent organisation ISCC. They guarantee that the amount of raw material used when manufacturing ECO+ medicine measures is supplied to the system in the form of fossil-free raw material.

By using independent mass balance certification we guarantee the mix of fossil and renewable raw materials will be followed up.



ECO+ wide medicine measure

30 ml medicine measure graduated each 5 ml.

The measures are made of polypropylene.

Accessorie: Lids .

Measure: staple/carton 50/2 000

Lid: bag/carton 500/2 000

Art no:

10300 nature

10301 yellow

10302 blue

10303 red

10305 green

10314 lid



Colour codes

Narrow and wide medicine measure 30 ml exist in different colours. These colours is part of a medication system where the following colours represent a dose occasion.

Transparent	= Breakfast
Yellow	= Lunch
Blue	= Dinner
Red	= Evening
Green	= Extra dose





Eco medicine measure

30 ml medicine measure graduated each ml.

The measures are made of polypropylene and the writable lids of white polyethylene.

Measure: staple/carton 90/9 000

Lid: bag/carton 500/9 000

Art no:

01350 nature

01351 yellow

01352 blue

01353 red

01354 green

00364 lid

Eco medicine measure - new package

30 ml medicine measure graduated each ml.

Measure: staple/carton 80/480

Art no:

01360 nature

01361 yellow

01362 blue

01363 red

01364 green



Medicine measure, narrow

30 ml medicine measure graduated each ml.

The measures are made of polypropylene and the writable lids of white polyethylene.

Measure: staple/carton 80/4 800

Lid: bag/carton 500/4 000

Art no:

01250 nature

01251 yellow

01252 blue

01253 red

01254 green

00264 lid

Medicine measure, narrow

30 ml medicine measure graduated each ml.

The measures are made of polypropylene and the writable lids of white polyethylene.

Measure: staple/carton 75/3 750

Lid: bag/carton 500/4 000

Art no:

J10823 nature

J10826 yellow

J10831 blue

J10828 red

J10833 green

00264 lid





Wide medicine measure

30 ml medicine measure graduated each 5 ml. The measures are made of polypropylene and the writable lids of white polyethylene.

Measure: staple/carton 50/2 000

Lid: bag/carton 500/2 000

Art no:

10300 nature

10301 yellow

10302 blue

10303 red

10305 green

10314 lid



Trays



Combination tray

Combination tray for 28 (7x4) medicine measures, narrow or wide. Made in polystyrene.

Size: 375 x 225 x 40 mm

15 pcs. per carton

Art no:

500



Tray for narrow measures

Medicine tray made of white polystyrene for 90 narrow measures.

Size: 495 x 377 mm.

3 pcs. per carton

Art no:

J10681



Dispenser - narrow Eco

Made of PMMA.

Holds 90 units.

Outside Dimensions:

395 x 50 mm

5 pcs. per carton.

Art no:

300650

Dispenser - narrow

Made of PMMA.

Holds 80 units.

Outside Dimensions:

395 x 50 mm

5 pcs. per carton.

Art no:

300450



Dispenser - wide

Made of PMMA.

Holds 50 units.

Outside Dimensions:

390 x 60 mm

5 pcs. per carton.

Art no:

300400



Dispenser for lids

Made of PMMA.

Outside Dimensions: 270 x 170 x 190 mm

5 pcs. per carton.

Art no:

300600



Tablet crushers and splitters



Tablet crusher, narrow

Crushes the tablets directly in the medicine measure.

10 pcs. per inner box

50 pcs. per carton

Art no:

300280



Tablet crusher, wide

Crushes the tablets directly in the medicine measure.

10 pcs. per inner box

50 pcs. per carton

Art no:

300320



Electric pill grinder Severo 3.0

Rechargeable and mobile. User-friendly and prevents RSI. Grinds in 8 seconds. Easy to clean. Silent running.

Weight: 1,5 kg

Dimensions: 80x70x217 mm

Art no:

300300

Severo Cup

Medicine cup to Severo pill grinder

Can be supplemented with rondels and writable lids.

Cup volume 30 ml

Size: staple/carton 105/2100 pcs

Lid: bag/carton 500/2000 pcs

Rondels: 2100 per box

Art no:

10306 cup

10314 lid

3000004 rondel



VitaCarry™ power pill grinder

Electric, rechargeable handheld pill grinder. Fast, easy and small enough to fit in your purse or pocket. The grinding cartridge is easily removed for safe cleaning. Rechargeable batteries and charger are included.

Size: 120 x 65 x 25 mm
12 pcs. per carton.

Art no:
300055



Tablet crusher

Tablet crusher made of polycarbonate. Crushes the tablets with a simple flick of the wrist.

8 pcs. per inner box
64 pcs. per carton

Art no:
300100



Pill puncher

Holds blister packs of medicine in place and easily removes the pill from hard-to-open packages.

Made of PP
Dimensions: 56 x 32 x 21 mm
10 pcs. per carton.

Art no:
300200



Tablet splitter

Tablet splitter made of polycarbonate. Splits the tablet, regardless of size, into two parts.

12 pcs. per inner box
120 pcs. per carton

Art no:
300350





Denture cup

Denture cup 270 ml, made of white polypropylene, with lid made of white polystyrene. Diameter/depth: 95/50 mm. 168 pcs. per carton.

Art no:

200270 cup w. lid



Hpm offer

Medical solutions in the areas of Surgery, Disposables and Custom.

Our mission is to develop, manufacture and sell innovative and effective solutions for the health care industry and thereby contribute to the continuous improvement of the same.



Version: 2018:2 2018-03-22

HpmTM

Hammarplast Medical AB

Headquarters

Box 2069, SE-351 02 Lidköping,
Sweden
Visiting address: Kartåsgatan 8
SE-531 02 Lidköping, Sweden

Branches

Hammarplast Medical AB, Bredaryd,
Västra Vägen 10
SE-330 10 Bredaryd, Sweden

Hammarplast Medical AS, Tallinn
Sõstramäe 8,
EE-114 15 Tallinn, Estonia

Tel. +46 (0)510 618 80

E-mail. info@hpm.se

Website www.hpm.se