

Surgery

Accessories for Tourniquets:
Cuffs and Roll-on cuffs



For better care!








Hpm[™]

Hpm™ Tourniquet cuffs, disposable



- ETO Sterilized
- 2 m long hose made of TPU
- Straight and Conical
- No Latex
- Wide range of sizes
- Color coded
- Customized stockinett included
- CPC male as standard
- Soft and easy to apply
- High quality material





Single disposable cuffs

Art. no.		Total length (cm)	Width (cm)	Total length (inch)	Width (inch)
15-110		30	8	12	3
15-120		38	10	15	4
15-130		46	10	18	4
15-140		61	10	24	4
15-150		76	13	30	5
15-160		86	13	34	5
15-180		109	13	43	5

Note! All sizes are valid with +/- 1 cm width and length

Hpm™ Tourniquet cuffs, disposable

Single disposable conical cuffs

Art. no.		Total length (cm)	Width (cm)	Total length (inch)	Width (inch)
75-140		61	10	24	4
75-150		76	10	30	4
75-160		86	11	34	4
75-170		109	13	43	5











Note! All sizes are valid with +/- 1 cm width and length



Tourniquet cuffs, autoclavable

- Many different sizes of both single and double cuffs.
- Both washable and autoclavable cuffs.
- No risk for latex allergy.
- Velcro and strap with plastic clamp for increased safety.

Single autoclavable cuffs









Art. no.		Total length (cm)	Width (cm)	Total length (inch)	Width (inch)
35-100		39	4	15,4	1,6
35-110		47	6	18,5	2,4
35-115		51	10	20,1	3,9
35-120		53	8	20,9	3,1
35-130		65	12	25,6	4,7
35-140		78	12	30,7	4,7
35-145		86	14	33,9	5,5
35-150		93	13	36,6	5,1
35-160		102	13	40,2	5,1
35-170		128	13	50,4	5,1

Note! All sizes are valid with +/- 1 cm width and length

Tourniquet cuffs, autoclavable



Double autoclavable cuffs

Art. no.		Total length (cm)	Width (cm)	Total length (inch)	Width (inch)
45-90		31	7	12,2	2,8
45-100		41	7	16,1	2,8
45-105		46	10	18,1	3,9
45-110		50	10	19,7	3,9
45-120		67	17	26,4	6,7
45-130		82	17	32,3	6,7
45-140		100	17	39,4	6,7
45-150		125	17	49,2	6,7

Note! All sizes are valid with +/- 1 cm width and length

Roll-on cuffs

Hpm™ roll-on cuff system

The Hpm™ roll-on cuff is a complete system in order to create and keep bloodless field in extremities during operations in bloodless field. The cuff is efficient, simple and safe to use and need no electricity or compressed air.

You can either use only the roll-on cuffs to empty the extremities of blood, or use them along with our usual cuffs to empty and maintain bloodless field.

The design secures an ideal pressure over the systolic pressure against the extremity within the working range of each cuff and can be sterilized.

Our roll-on cuffs are made of silicone thus no risk for latex allergy.


Working areas for the cuffs:

- Orthopaedic surgery
- Vein surgery
- Leg and arm surgery
- Emergency care



Roll-on cuffs

Art. no.	Patient limb size
2531	25-31 cm
3038	30-38 cm
3746	37-46 cm
4556	45-56 cm
5565	55-65 cm
6480	64-80 cm



Roll-on cuffs



Vacuum pump
Art. no. 9007



Handpump (bladder)
including manometer
Art. no. 9005



Special wedge
Rubber, no risk for latex allergy
Art. no. 9001



Measuring tape
Art. no. 9004



Hose with needle and nipple
Acc. to pump art.no. 9005
Art. no. 9003



Roll-on cuff set Standard

The set consists of:
2531, 3038, 3746, 4556, 5565,
2 x 9001, 9003, 9004, 9005 and 9007.

Art. no. 9100

Roll-on cuff set Leg

The set consists of:
3038, 3746, 4556, 5565, 6480, 2 x 9001,
9003, 9004, 9005 and 9007.

Art. no. 9200

Hpm offer

Medical solutions in the areas of Surgery, Disposables and Custom.

Our mission is to develop, manufacture and sell innovative and effective solutions for the health care industry and thereby contribute to the continuous improvement of the same.



Version: 2020-03-24 WWW

HpmTM

Hammarplast Medical AB

Headquarters

Box 2069, SE-351 02 Lidköping,
Sweden
Visiting address: Kartåsgatan 8
SE-531 02 Lidköping, Sweden

Branches

Hammarplast Medical AB, Bredaryd,
Västra Vägen 10
SE-330 10 Bredaryd, Sweden

Hammarplast Medical AS, Tallinn
Sõstramäe 8,
EE-114 15 Tallinn, Estonia

Tel. +46 (0)510 618 80

E-mail. info@hpm.se

Website www.hpm.se