

Ear Speculum

NEW!



Hpm[™]

For better care!

Ear speculum disposable



Hpm™ disposable ear speculums are hygienic and safe.

- Available in two different sizes: ø 2.5 mm and ø 4 mm.
- Disposable ear speculums are hygienic and prevents from cross-contamination.
- Secure fit
- No sharp edges
- Available in colour black
- Made of recycleable polypropylene, do not contain BPA, phtaltes, latex or PVC.
- Hpm™ ear speculums are CE-marked and classified as class 1.
- For single use

REF nr	Size mm	Colour	Package	EAN code
31110	ø 2,5	black	250 pcs/box	7340077904058
31120	ø 4	black	250 pcs/box	7340077904072



Hpm™ ear speculum fits the most commonly used otoscopes.

Otoscope	
Heine Alpha+	KaWe Combilight F.O
Heine K180	KaWe Eurolight F.O
Heine BETA 400, 800	Riester Ri-mini
Heine MINI 2000 + FO	Riester Ri Scope L1
Heine MINI 3000+FO+ 3000	Riester Ri Scope L2
Heine BETA 200	Riester E-scope
KaWe Piccolight C	



Hpm™
For better care!

Hpm offer

Medical solutions in the areas of
Surgery, Disposables and Custom.

Our mission is to develop, manufacture
and sell innovative and effective
solutions for the health care industry
and thereby contribute to the
continuous improvement of the same.



Version: Ear speculum_ 2019:2 2019-03-28

HpmTM

Hammarplast Medical AB

Headquarters

Box 2069, SE-351 02 Lidköping,
Sweden
Visiting address: Kartåsgatan 8
SE-531 02 Lidköping, Sweden

Branches

Hammarplast Medical AB, Bredaryd,
Västra Vägen 10
SE-330 10 Bredaryd, Sweden

Hammarplast Medical AS, Tallinn
Sõstramäe 8,
EE-114 15 Tallinn, Estonia

Tel. +46 (0)510 618 80

E-mail. info@hpm.se

Website www.hpm.se