


- ISO 13485
- ISO 9001
- ISO 14001

## Medicine measure, wide 30 ml (2000)



## Product information

Ref.nr	Name	Colour	Diameter	Height	Volume
10300	Medicine measure, wide	Transparent	ø 47mm	32 mm	30 ml
10301	Medicine measure, wide	Yellow	ø 47mm	32 mm	30 ml
10302	Medicine measure, wide	Blue	ø 47mm	32 mm	30 ml
10303	Medicine measure, wide	Red	ø 47mm	32 mm	30 ml
10305	Medicine measure, wide	Green	ø 47mm	32 mm	30 ml


## Packaging

Ref nr	Inner package (sleeve)	Outer Package, case	Cases per pallet	Measurement Case	Package pallet/height	Weight package	Weight pallet	Barcode
10300	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077901217
10301	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077901224
10302	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077901231
10303	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077901248
10305	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077901262

All our outer packages are marked with FSC logotype. The FSC-logo are a trusted sign of responsible use of the planet's forest resources.


# Information

- Disposable product for easy and quick handling, increased hygien and reduced risk of contamination.
- Medicine measure suitable for dispensing both liquid and dry medications (tablets)
- Graduated every 5 milliliter up to 30 milliliter. We believe an exact scale is very important when you distribute medicine to a patient.
- High quality Medicine measure made of recycleable PP.
- Classified as a class 1M product. Class 1M means that there is a measurement function with a scale on the product. All 1M products must have Notified Body mark and CE mark on the product. Our Notified Body number 0413 implies that our Notified Body, Intertek, inspects our calibration system of the medicine cups and controls that we follow the manufacturing rules of the product according to the Medical Device regulations.
- Made of polypropylene
- Accessories:
  - Writeble lids made of PE, ref nr 10314
  - Tablet crusher ref.nr. 300320
  - Dispenser ref.nr. 300400



# Labels

Coloured line for faster identification of right product.





Article name, colour and pcs in package

Symbols

EAN code

# Symbols



Product conforms with applicable EC directives. Notified body number, Inspected yearly by Intertek.



Packaging from responsible resources