



- ISO 13485
- ISO 9001
- ISO 14001

ECO Medicine measure, narrow 30 ml 9000-p







Product information

Ref.nr	Name	Colour	Diameter	Height	Volume			
01350	ECO Medicine measure, narrow	Transparent	ø 38mm	42 mm	30 ml			
01351	ECO Medicine measure, narrow	Yellow	ø 38mm	42 mm	30 ml			
01352	ECO Medicine measure, narrow	Blue	ø 38mm	42 mm	30 ml			
01353	ECO Medicine measure, narrow	Red	ø 38mm	42 mm	30 ml			
01354	ECO Medicine measure, narrow	Green	ø 38mm	42 mm	30 ml			

Packaging


Ref nr	Inner package (sleeve)	Outer Package, case	Cartons per pallet	Size package	Package pallet/height	Weight package	Weight pallet	Barcode
01350	90	9000	30	390x390x415 mm	228 cm	9,5 kg	307 kg	7340077903532
01351	90	9000	30	390x390x415 mm	228 cm	9,5 kg	307 kg	7340077903488
01352	90	9000	30	390x390x415 mm	228 cm	9,5 kg	307 kg	7340077903495
01353	90	9000	30	390x390x415 mm	228 cm	9,5 kg	307 kg	7340077903556
01354	90	9000	30	390x390x415 mm	228 cm	9,5 kg	307 kg	7340077903549

All our outer packages are marked with FSC logotype. The FSC-logo are a trusted sign of responsible use of the planet's forest resources.



Information


- Disposable product for easy and quick handling, increased hygien and reduced risk of contamination.
- Medicine measure suitable for dispensing both liquid and dry medications (tablets)
- Graduated each milliliter up to 30 millilter. We belive an exact scale is very important when you distribute medicine to a patient. Note that the medicine description always prescribes an exact dose. If the product is not marked with Notified Body number the scale can not be trusted.
- Classified as a class Im product. Class Im means that there is a measurement function with a scale on the product. All Im products must have Notified Body mark and CE mark on the product. The Notified Body Number (N.B) signifies that the product is following the regulations of class Im, The Notified Body is inspecting the manufacturing company yearly regarding the calibration of the scale
- High quality medicine measure made of PP.
- Accesories:
 - Writeble lids made of PE, ref nr 00364
 - Tablet crusher ref nr 300280
 - Dispenser Ref.nr 300650



Labels



Symbols



Packaging from responsible resources

Hammarplast Medical AB Kartåsgatan 8 53440 Lidköping

info@hpm.se www.hpm.se

For better care!