

ECO+ Medicine measure, wide 30 ml (2000)



Product information

Ref.nr	Name	Colour	Diameter	Height	Volume
02650	ECO+ Medicine measure, wide	Transparent	ø 47mm	32 mm	30 ml
02651	ECO+ Medicine measure, wide	Yellow	ø 47mm	32 mm	30 ml
02652	ECO+ Medicine measure, wide	Blue	ø 47mm	32 mm	30 ml
02653	ECO+ Medicine measure, wide	Red	ø 47mm	32 mm	30 ml
02654	ECO+ Medicine measure, wide	Green	ø 47mm	32 mm	30 ml


Packaging

Ref nr	Inner package (sleeve)	Outer Package, case	Cases per pallet	Size package	Package pallet/height	Weight package	Weight pallet	Barcode
02650	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077904003
02651	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077904010
02652	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077904027
02653	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077904034
02654	50	2000	56	370x234x375 mm	228 cm	4,5 kg	274 kg	7340077904041

All our outer packages are marked with FSC logotype. The FSC-logo are a trusted sign of responsible use of the planet's forest resources.


Information

- Disposable product for easy and quick handling, increased hygien and reduced risk of contamination.
- Medicine measure suitable for dispensing both liquid and dry medications (tablets)
- Graduated every 5 milliliter up to 30 milliliter. We believe an exact scale is very important when you distribute medicine to a patient.
- High quality Medicine measure made of recycleable PP.
- Classified as a class 1M product. Class 1M means that there is a measurement function with a scale on the product. All 1M products must have Notified Body mark and CE mark on the product. Our Notified Body number 0413 implies that our Notified Body, Intertek, inspects our calibration system of the medicine cups and controls that we follow the manufacturing rules of the product according to the Medical Device regulations.
- Made of polypropylene
- Accessories:
 - Writeble lids made of PE, ref nr 10314
 - Tablet crusher ref.nr. 300320
 - Dispenser ref.nr. 300400



Labels

Coloured line for faster identification of right product.





Article name, colour and pcs in package

Symbols

EAN code

Symbols



Product conforms with applicable EC directives. Notified body number, Inspected yearly by Intertek.



Made of renewable raw material, inspected yearly by the independent organisation ISCC.



Packaging from responsible resources

Hammarplast Medical AB
Kartåsgatan 8
53440 Lidköping

info@hpm.se
www.hpm.se

For better care!