

Medicine measure narrow 30 ml



Product information

Ref.nr	Name	Colour	Diameter	Height	Volume
01250	Medicine measure, narrow	Transparent	ø 38mm	42 mm	30 ml
01251	Medicine measure, narrow	Yellow	ø 38mm	42 mm	30 ml
01253	Medicine measure, narrow	Blue	ø 38mm	42 mm	30 ml
01254	Medicine measure, narrow	Red	ø 38mm	42 mm	30 ml
01255	Medicine measure, narrow	Green	ø 38mm	42 mm	30 ml

Packaging

Ref nr	Inner package (sleeve) pcs	Outer Package, case pcs	Cases per pallet	Size package	Package pallet/height	Weight package	Weight pallet	Barcode
01250	80	4800	56	370 x 234 x 375 mm	232cm	8 kg	470 kg	7340077900036
01251	80	4800	56	370 x 234 x 375 mm	232 cm	8 kg	470 kg	7340077900043
01252	80	4800	56	370 x 234 x 375 mm	232 cm	8 kg	470 kg	7340077900050
01253	80	4800	56	370 x 234 x 375 mm	232 cm	8 kg	470 kg	7340077900067
01254	80	4800	56	370 x 234 x 375 mm	232 cm	8 kg	470 kg	7340077900074

All our outer packages are marked with FSC logotype. The FSC-logo are a trusted sign of responsible use of the planet's forest resources.