


- ISO 13485
- ISO 9001
- ISO 14001

Medicine measure narrow 30 ml



Product information

Ref.nr	Name	Colour	Diameter	Height	Volume
01250	Medicine measure, narrow	Transparent	ø 38mm	42 mm	30 ml
01251	Medicine measure, narrow	Yellow	ø 38mm	42 mm	30 ml
01253	Medicine measure, narrow	Blue	ø 38mm	42 mm	30 ml
01254	Medicine measure, narrow	Red	ø 38mm	42 mm	30 ml
01255	Medicine measure, narrow	Green	ø 38mm	42 mm	30 ml


Packaging

Ref nr	Inner package (sleeve) pcs	Outer Package, case pcs	Cases per pallet	Size package	Package pallet/height	Weight package	Weight pallet	Barcode
01250	80	4800	56	378x241x382 mm	232cm	8 kg	470 kg	7340077900036
01251	80	4800	56	378x241x382 mm	232 cm	8 kg	470 kg	7340077900043
01252	80	4800	56	378x241x382 mm	232 cm	8 kg	470 kg	7340077900050
01253	80	4800	56	378x241x382 mm	232 cm	8 kg	470 kg	7340077900067
01254	80	4800	56	378x241x382 mm	232 cm	8 kg	470 kg	7340077900074

All our outer packages are marked with FSC logotype. The FSC-logo are a trusted sign of responsible use of the planet's forest resources.


Information

- Disposable product for easy and quick handling, increased hygien and reduced risk of contamination.
- Medicine measure suitable for dispensing both liquid and dry medications (tablets)
- Graduated each milliliter up to 30 milliliter. We believe an exact scale is very important when you distribute medicine to a patient.
- High quality Medicine measure made of recyclable PP.
- Classified as a class 1M product. Class 1M means that there is a measurement function with a scale on the product. All 1M products must have Notified Body mark and CE mark on the product. Our Notified Body number 0413 implies that our Notified Body, Intertek, inspects our calibration system of the medicine cups and controls that we follow the manufacturing rules of the product according to the Medical Device regulations.
- Made of polypropylene
- Accessories:
 - Writeble lids made of PE, ref nr 00264
 - Tablet crusher ref.nr 300280
 - Dispenser ref.nr 300450



Labels

Symbols



Article name, colour and pcs in package

EAN code

Symbols

	Manufacturer		Do not re-use		Product conforms with applicable EC directives. Notified body number, Inspected yearly by Intertek.		Packaging from responsible resources
	Keep dry		Referens number		Made of renewable raw material, inspected yearly by the independent organisation ISCC.		
	Made of recyclable material		Lot number				

Hammarplast Medical AB
Kartåsgatan 8
53440 Lidköping

info@hpm.se
www.hpm.se

For better care!