New!


Hpm[™] medicine measure, made of renewable raw material.

The Hpm[™] vision calls for us to support the transformation towards a more sustainable society. One of our main goals is to increase production and sales of products made of renewable materials in order to reduce our carbon footprint.

Launching ECO+ medicine measure is one of our important contributions for a sustainable world.







Product information medicine measure 9000 pcs

- Made of renewable raw material.
- Disposable product for esy and quick handling, increased hygien and reduced risk of contamination.
- Medicine measure suitable for dispensing both liquid and dry medications (tablets)
- Graduated each milliliter up to 30 millilter. The scale is clear and easy to read for measuring.
- High quality medicine measure.
- Classified as a class 1M product. Class 1M means that there is a measurement function with a scale on the product. All 1M products must have Notified Body mark and CE mark on the product. Our Notofied Body number 0413 implies that our Notified Body, Intertek, inspects our calibration system of the medicine cups and controls that we follow the manufacturing rules of the product according to the Medical Device regulations.
- Manufactured in Sweden
- Accessories:
 - Writeble lids, ref nr 00364

INFORMATION	ART.NO.	COLOR	BARCODE
Packaging:	02350	Transparent	7340077903907
ECO+ medicine measure, narrow:	02351	Yellow	7340077903914
9000 pcs. per outer	02352	Blue	7340077903921
100 sleeves per inner	02353	Red	7340077903938
90 medicine measures per sleeve	02354	Green	7340077903945

