


New! Hpm™ ear speculum , made of renewable raw material.

The Hpm™ vision calls for us to support the transformation towards a more sustainable society. One of our main goals is to increase production and sales of products made of renewable materials in order to reduce our carbon footprint.

Launching ECO+ ear speculum is one of our important contributions for a sustainable world.



ECO+ Ear Speculum



Hpm™
For better care!

ECO+ Ear speculum



Hpm™ disposable ECO+ ear speculums are hygienic and safe

- Hpm™ ECO+ ear speculums are made of renewable raw material. Hpm™ is certified by the independent organisation ISCC and will be audited yearly.
- Available in two different sizes: \varnothing 2.5 mm and \varnothing 4 mm.
- Disposable ear speculums are hygienic and prevents from cross-contamination.
- Secure fit
- No sharp edges
- Available in colour black
- Made of recycleable polypropylene, do not contain BPA, phtaltes, latex or PVC.
- Hpm™ ear speculums are CE-marked and classified as

REF nr	Size mm	Colour	Package	EAN code
31210	ø 2,5	black	250 pcs/box	7340077904096
31220	ø 4	black	250 pcs/box	7340077904102



Hpm™ ear speculum fits the most commonly used otoscopes.

Otoscope

Heine Alpha+

Heine K180

Heine BETA 400, 800

Heine MINI 2000 + FO

Heine MINI 3000+FO+ 3000

Heine BETA 200

KaWe Piccolight C

KaWe Combilight F.O

KaWe Eurolight F.O

Riester Ri-mini

Riester Ri Scope L1

Riester Ri Scope L2

Riester E-scope

Hpm™ care even more about the environment,

FSC-marked carton

By only using FSC-labeled outer cartons we take additional responsibility for our environment. The FSC mark guarantees that the product is made of forest that is FSC-certified and responsible.

Hammarplast Medical AB
Box 2069, 531 02 Lidköping
+46 (0)510 618 80
hpm.se info@hpm.se

Hpm™



For better care!