

Medicine Measure, narrow blue

The wide medicine measure 30 ml are a class 1 M product and is manufactured according to the Medical Device regulation 93/42/EEC and its purpose is for dispensing both oral medication and tablets. The product is tested to fulfill all demands according to regulation. The narrow medicine cup can be used in medicine trays or together with a tablet crusher.

Our narrow medicine measure 30ml are a part of a medicine distributing system. The colours represent different dispensing occasions as follows:

Transparent	= Breakfast
Yellow	= Lunch
Blue	= Dinner
Red	= Evening
Green	= Extra dose

You can get lids to all narrow medicine measures so that you can protect the medication and you can also write on the lids.

Diameter:	38 mm
Height:	42 mm
Volume:	30 ml
Color:	Blue
Item no:	01252
Packaging:	stack/carton: 80/4 800



For better care!

HpmTM