



Eco medicine measure

the new improved product

Eco medicine measures are designed to meet the patients' and nurses' needs:

- 30 ml medicine measure graduated each ml. The scale is clear, reliable and easy to read for measuring.
- For your patients' safety our Eco medicine measures are CE-marked and Notified Body marked according to Medical Device Directive, MDD 93/42/ECC.
- Made in transparent polypropylene PP, BPA free.
- More comfortable drinking edge, smooth to the lips.
- Improved grip function for your patients.
- Stable top edge for easy handling.
- Available in 5 colors representing different dispensing occasions.
- Manufactured in Sweden.
- Can be completed with a writable lid.

INFORMATION

Packaging :

Eco medicine measure:

9000 pcs. per outer

100 sleeves per inner

90 medicine measures per sleeve

Lids:

9000 per outer

ART.NO.

FÄRG

EAN

01350

Transparent

7340077903532

01351

Yellow

7340077903488

01352

Blue

7340077903495

01353

Red

7340077903556

01354

Green

7340077903549

00364

White lid

7340077903525

Disposable product



CE marked

The Notified Body Number signifies that the product is following the regulations of class 1M, where the Notified Body is inspecting the manufacturing company regarding the calibration of the scale.



For better care!

Hpm™

Hammarplast Medical AB
Box 2069, SE-531 02
Lidköping, Sweden
Tel: +46 (0)510 618 80
info@hpm.se
www.hpm.se