

New!

Hpm™ medicine measure, made of renewable raw material.

The Hpm™ vision calls for us to support the transformation towards a more sustainable society. One of our main goals is to increase production and sales of products made of renewable materials in order to reduce our carbon footprint.

Launching ECO+ medicine measure is one of our important contributions for a sustainable world.



ECO+ Medicine measure



HpmTM
For better care!



Medicine measure made of renewable raw material

Product information medicine measure wide

- Made of renewable raw material.
- Disposable product for easy and quick handling, increased hygien and reduced risk of contamination.
- Medicine measure suitable for dispensing both liquid and dry medications (tablets)
- Graduated every 5 milliliter up to 30 milliliter. The scale is clear and easy to read for measuring.
- High quality medicine measure.
- Classified as a class 1M product. Class 1M means that there is a measurement function with a scale on the product. All 1M products must have Notified Body mark and CE mark on the product. Our Notified Body number 0413 implies that our Notified Body, Intertek, inspects our calibration system of the medicine cups and controls that we follow the manufacturing rules of the product according to the Medical Device regulations.
- Manufactured in Sweden
- Accessories:
 - Writeble lids, ref nr 10314
 - Tablet crusher, narrow measures ref.nr. 300320

INFORMATION

Packaging:

ECO+ medicine measure:

2000 pcs. per outer

40 sleeves per inner

50 medicine measures per sleeve

ART.NO.

02650

02651

02652

02653

02654

COLOR

Transparent

Yellow

Blue

Red

Green

BARCODE

7340077904003

7340077904010

7340077904027

7340077904034

7340077904041

HpmTM
For better care!

