## ARTNO. 300400 (wide), 300450 (narrow)

## Dispenser for cups

## **INFORMATION**

Transparent dispenser for wide and narrow medicine cups.

Organize your medicine cups with a dispenser that can be mounted on a wall or on a medicine trolley. Removable cap on the top of the dispenser allows for convenient loading and keeps dust and dirt away. Transparent allows an inventory check at a glance. The CE marked dispenser is delivered with two screws for easy mounting. Made in PMMA.

Wide, 300400, holds 50 units, Narrow, 300450, holds 80 units.

## SIZE (w x h)

Wide, 300400 60 x 390 mm

Narrow, 300450 50 x 390 mm





